

Optimal 5zon®

Flexible mattress replacement for pressure ulcer prevention







Optimal 5zon® with cover Olivia/Evac

FACTS	OPTIONS
 Core: HR-foam 65+50 kg/m3 Height: 14 cm Function: as an aid for prevention and treatment up to and including category 2⁽¹⁾ Recommended user weight: up to 180 kg Fire: EN 597-1, -2 Warranty: 3 years Cleaning of cover: wipe with cleaning agent and/or disinfectants. Machine wash max 95 °C, tumble drying 	• Inner cover SIZES 80/85/90/105/120x200 cm and 80/85/90x210 cm. If you need any other sizes please contact your distributor.
and/or disinfectants.	cm. If you need any other sizes please

Optimal 5zon®

Optimal 5zon® is intended for all healthcare environments as an aid in the prevention and treatment of pressure ulcers up to and including category 2⁽¹⁾. It is a flexible mattress replacement that is placed directly on the bed and replaces both the standard mattress and mattress pad.

Optimal 5zon® has good pressure redistribution properties⁽²⁾ and is also equipped with our patented⁽⁵⁾ DuoCore® sit function. The sit function has been designed to provide support when seated, while also offering good comfort in lying position.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene, quality and safety⁽³⁾⁽⁶⁾. All covers are manufactured from modern, recyclable materials and are vapour permeable⁽⁴⁾, which contributes to reduce the risk of skin maceration.

REFERENCES

- (1) European Pressure Ulcer Advisory Panel, National Pressure Injury Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers/Injuries: Clinical Practice Guideline. The International Guideline. Emily Haesler (Ed). EPUAP/NPIAP/PPPIA: 2019
- (2) Examination of the shear force and the pressure relief capability, Berlincert 2009.
- (3) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019.
- (4) SS 876 00 20:2017, Sjukvårdstextil-Madrasser-Specifikationer och krav.
- (5) Patentnummer 008522741.
- (6) SS 876 00 20:2017, Healthcare textiles Mattresses Specifications and requirements













Always read the instructions for use prior to usage of the product.

Care of Sweden, Optimal 5zon and DuoCore are trademarks of Care of Sweden AB. © Care of Sweden AB, 2021.

All rights reserved. Last update: 2021-06-03.

FEATURES & CHARACTERISTICS

Pressure

Elastic material and a zone and module system gives pressure equalisation.

Shear forces

The mattress is specially designed to reduce the effects of shear forces.

Sit function

Our patented⁽⁵⁾ sit function DuoCore[®] is designed to give support when sitting and also to provide good comfort when lying.

Moisture

The hygienic cover is vapour permeable⁽⁴⁾, which contributes to reduce the risk of skin maceration.

Hygiene

The cover is easy to clean and has a good hygienic standard⁽³⁾⁽⁶⁾.

AVAILABLE COVERS

Hugiene cover Olivia

Drip shield zip on three sides.

Hygiene cover Olivia H

Drip shield zip on three sides and with handles for transport and moving.

Hygiene cover Nordic

Welded seams, drip shield zip on three sides and handles on one side.

Hygiene cover Olivia/Evac

A 2-piece cover with handles for evacuation. It has a drip shield zipper on all sides.